



400.50 Commode shower char with rear wheels Ø 600 mm



Categoria Prodotto: [IHS](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/400-50-commode-shower-char-with-rear-wheels-600-mm/>

Descrizione Prodotto

Structure made of anodized aluminium tube powder varnished. Ring shaped removable padded PU high density seat with antiscratch finishing. Padded PU armrests adjustable in height and removable. Thin padded backrest adjustable in height footplates removable for easy access of the user, revolving and adjustable in flex extension. Two wheels from Ø 125 mm and two rear wheels from Ø 600 mm with nylon support and parking brakes. Handle for the attendant. Plastic, backwards removable evacuation pan. Demountable frame for transport. All the screws and fixing components are made of Inox Steel 304 rustproof. All the materials used are hypoallergenic and atoxic.

Load capacity kg 120.

Requirements following Rule EN ISO 12183.

Rear wheels Ø 600 mm