



500.10 FOLDING LIGHT WHEELCHAIR



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/500-10-folding-light-wheelchair/>

Descrizione Prodotto

Structure made of varnished aluminium, double cross, padded seat, backrest and armrests, calfrest strap adjustable and removable, upsetting and removable armrests with side panels, seat adjustable in depth of 3 cm.

Footrests with angle adjustment turning outwards and inwards, removable by quick

release. Rear pneumatic wheels removable by quick release with parking brakes and handrim for selfpropelling. Bag on the rear of the backrest.

Front driving wheels Ø 160 puncture proof.

Available seating sizes: cm 41-46 Weight kg 15. Load capacity kg 120.

Requirements following Rule EN ISO 12183. EN ISO 12183.