



700.10 FOLDING STANDARD WHEELCHAIR



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/700-10-folding-standard-wheelchair/>

Descrizione Prodotto

Structure made of varnished steel tube, double cross, padded seat, backrest and armrests. Armrests backwards upsetting and folding until laying on the bottom of the backrest, calfrest strap adjustable and removable.

Footrests revolving, removable and adjustable. Handles for the attendant, rear

stirrups to help the step climbing, rear pneus wheels with handrim for selfpropelling and parking brakes.

Front wheels Ø 200x32 puncture proof, bag on the rear of the backrest.

Available seating sizes: cm 41-46. Weight kg 18. Load capacity kg 120.

Requirements following Rule EN ISO 12183.