



## 10.70XXL NURSING 300



**Categoria Prodotto:** [Extra Large Solutions](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/nursing-300/>

### Descrizione Prodotto

**Bed with 4 sections.** The frame is made of powder varnished steel shapes, electrowelded and reinforced bottom with 4 articulated sections with double actuators, the backrest is with rear transfer to maintain the user in a stable position

during the inclination of the backrest in order to favour a reciprocal inclination with the mattress inserted, reclining legs section until 45° with feet area to follow on side pantographs.

Feeding of the movements by linear actuators at 24V with handset with push buttons on spiral cable, buffer batteries inserted in the power module for emergency in case of lack of power supply, 4 legs at the base setted up with holes adjustable from cm 40 to 52, Nursing head and feet panels in MDF covered on the 2 sides in colour walnut and sand.

**Dim. cm 203x120x40+12. Max. load kg. 300**