



VARIO PLUS 12.17E



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/vario-plus-12-17e/>

Descrizione Prodotto

NEW VARIO PLUS ELECTRIC MULTIFUNCTION

Multifunction wheelchair with functions electroassisted feeded by battery, functions

on handset. Electrically tilting seating, until 45°, and reverse tilting until 10°, and handly adjustable in height until + 6 cm and in depth until + 12 cm.

Electrical independent reclining legrests and removable; opening footplates handly adjustable in flexion and extension and padded adjustable calfrests. Removable side panels, adjustable in width and depth. Padded reclining armrests, adjustable in height and inclination and equipped with system which allows to follow the inclination movement of the backrest during this function.

Padded shaped seat applied on rigid bottom and removable with velcro. Electric reclining padded wrapping shaped backrest on rigid aluminium bottom shaped and removable with velcro; padded contoured headrest adjustable in height and depth on rotating base and adjustable, rotating attendant handles with levers to operate the wheel brakes.

Rear pneumatic wheels Ø 300 mm with brakes, removable and revolving antitipping castors, front pivoting steering wheels Ø 200 mm. Max. load capacity 120 Kg.

* Available also with rear puncture proof wheels Ø 600 mm as detailed.