



VARIO PLUS 12.17T



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/vario-plus-12-17t-2/>

Descrizione Prodotto

WHEELCHAIR WITH TILTING FRAME SERVO-ASSISTED

Suitable for users with limited mobility which need continuous changements of the

position to compensate the breathy difficulties, post surgery, circulatory, articular and/or of physical strenght imposed by dystonia and/or spasms. Seating adjustable in depth by reciprocal articulations.

Armrests removable and adjustable in height. Footplates reclining, adjustable and removable. Tilting until 45° by gas pump, backrest inclination until 60° by gas pump. Padded seat fixed on the bottom with fabric removable on velcro. Padded wrapping backrest adjustable by slings. Padded headrest semiinvasive adjustable in height and depth on rotating base.

Available in two seating width cm 40 - 45. Max load 120 kg.